



DATA SHEET

CanTab®

pH at 20% Solution	3.5–5.5
Conductivity, micomhos	<100
Dextrose Equivalent, %	95.0 – 99.9
Light Transmittance, %	80 – 100
Moisture, 100%	7.5 – 9.0
Sieve Analysis on a 20 Mesh	0 – 3
Sieve Analysis Thru a 100 Mesh	0 – 25
Total Count/ gm Composite	Less than 1000
Mold Count / gm Composite	Less than 50
Yeast Count / gm Composite	Less than 500

PACKAGING AND STORAGE

Net weight 50 lbs. in 3-ply paper bag, with polyethylene lining. Keep product away from water, humidity and high temperature to avoid moisture absorption